

## SQUARE WAVE CATHODIC STRIPPING VOLTAMMETRY ADSORPTIVE FOR NICKEL AND COBALT ANALYSIS

Saryati\*, Rukihati and Sumardjo

Technology Center for Nuclear Industrial Material - National Nuclear Energy Agency (BATAN), Jakarta

Received 9 March 2006; Accepted 30 May 2006

### ABSTRACT

The adsorptive stripping voltammetric determination of Ni and Co based on adsorption of the Ni/Co and dimethylglyoxime (DMG) complex on a hanging mercury drop electrode is studied. The reduction current of the adsorbed DMG complex is measured by square wave cathodic stripping voltammetry method. The effect of various parameters such as ligand concentration, pH of supporting electrolytic, adsorption potential and adsorption time on the current peak of Ni and Co voltammogram were studied. Optimum condition of this method are supporting electrolyte pH 9, DMG concentration  $5 \times 10^{-4}$  M, adsorption potential -0.7 V vs Ag/AgCl and adsorption time 180 second. A linear relationship between the current peak and Ni or Co concentration was obtained in the range 5 – 30 ng/mL and the detection limit 0.6 ng/ml for both Ni and Co. The recovery of Ni and Co were 98.11-104.17% using standard biological materials with RSD 2.59 – 10.37%. Based on "t" test can be conclude that the result are nearly equal to the standard reference material.

**Keywords:** adsorptive stripping voltammetric, dimethylglyoxime complex, nickel, cobalt.

### INTRODUCTION

In a row with the awareness increase of life, environment and materials quality, the need of quality analysis method are rise. Industrial development caused heavy metal pollution. So that heavy metal analysis is one interesting study in life and environment science

Heavy metals are a great significance in ecochemistry and ecotoxicology because of their toxicity at low levels and tendency to accumulate in human organs. The major intake of toxic metals by humans is through food. Vegetable foods are relative more exposed to environmental pollution than cereals. Leafy vegetables are most exposed to environmental pollution because of their large surface area. In many countries vegetables are a major part of diet and its important to monitor their heavy metal content. Cobalt and nickel are very important heavy metal in biochemical and psychological process, but accumulate and deficiency cause health interference [1]

There is many instrument and method that used in heavy metal analysis, voltammetry method is one of it. The advantage of voltammetry method is a major element not interference in the minor element determination, so that the sample preparation is faster, because it doesn't need element separation. In other hand voltammetry method can be also used to determine ion in wide concentration range, from ng/ML(ppb) until  $\mu\text{g}/\text{ML}$  (ppm) level. For height concentration ( $\geq$  ppm) direct analysis by using differential pulse voltammetry (DPV) or *square wave voltammetry* (SWV) is possible. But for low concentration the DPV and SWV technique can't be used except it combining with accumulation step, that know as

stripping voltammetry, anodic stripping voltammetry, cathodic stripping voltammetry or adsorptive stripping voltammetry.

In a media that can't form complex compound Ni(II)/Ni(I) or Co(II)/Co(I) pairs electrode reaction are very irreversible, so its hard to determinate using anodic/cathodic stripping voltammetry, because of that Ni and Co determination was done by using adsorptive stripping voltammetry with complexing agent. The complexing agent increase the sensitivity method. Some researcher used dimethylglyoxime (DMG) as complexing agent [3-5]. Neiman [3] has studied Ni adsorptive stripping voltammetry in tetraethyl amine and ammonium chloride media with Ep Co -1,04 V/SCE. Vos *et al* [4] has been determined Ni and Co in rain water in the ammonium buffer media with Ep Ni - 1.04 V/SCE and Ep Co -1.195 V/SCE. Thunus [5] has determined Ni and Co in blood with pH 9 ammonium buffer media. The several literature used DPV in cathodic current from Ni and Co -DMG complex reduction. The optimum condition of this method commonly determinate by experimentally with take a note of sensitivity, reproductivity and peak voltammogram form.

For reliable analysis by voltammetry, the analyte must be in the ionic state. In voltammetric analysis where usually associated with the incomplete destruction of organic matter during digestion. Wet digestion in the mixture nitric acid perchloric acid and hydrogen peroxide that follow by heating, has proven successfully to determinate metal ion in biological sample [2-5].

This paper reported the optimization of Ni and Co determination simultaneously with adsorptive cathodic

\* Corresponding author.

square wave stripping voltammetric method by using DMG as complexing agent and application of this method to determine Ni and Co in Apple leave SRM 1515, Bovine liver SRM 1577B, Oyster Tissue SRM 1566b, Sargasso CRM no 9, Copepod RMMAA-1 and local vegetables

## EXPERIMENTAL SECTION


### Instrumentation and reagents

The voltammograms were recorded on a EG&G PAR model 384B polarographic analyzer. The working electrode was a EG&G PAR model 303 A static mercury drop electrode and hanging mercury drop electrode (HMDE) was employed. Ag/AgCl (saturated KCL) and platinum wire as reference and auxiliary electrode, respectively. A XY recorder for Houston instrument was used.

All reagent were of analytical-reagent grade, except a stock solution of 10000 µg/mL Ni and Co from MERCK was used. All solution were prepared from distilled water.

### Procedure

Voltammetric stripping experiments were described in detail below. The supporting electrolyte (10 ml) containing 0.1 M ammonium tartaric buffer,  $5 \cdot 10^{-4}$  M DMG transferred into the voltammetric cell and purged with nitrogen gas for 4 min. The pre-concentration potential at  $-0.7$  V was applied to a fresh mercury drop while the solution was stirred. Stirring was stopped and then after 5 s voltammogram recorded by applying a negative scan (at 200 mV/s) from  $-0.7$  to 1.3 V using the square wave mode. After the ground voltammogram has been obtained, the adsorptive stripping experiment was repeated with a new drop for the addition of sample. Nitrogen was then purge through the solution for about 1 min. Quantification was done by standard addition using three spikes for each element.




**Fig 1.** The effect of DMG concentration on current peak (Ip) Ni and Co. Ni concentration = Co concentration = 20 ng/mL. Deposition time 120 s in tartaric buffer pH 9 media

A covered quartz cell containing a 0.4 g sample powder and 2 mL concentrated  $\text{HNO}_3$ , 0.5 mL perchloric acid were added and placed on the hot plate and keep slightly boiling until dissolution of the sample was complete. If dissolution was in complete, the mixture was cooled to  $80$  °C and 1 mL concentrated  $\text{HNO}_3$ , 0.4 mL  $\text{H}_2\text{O}_2$  30% was added. The mixture was heated to dryness to yield a whitish residue. The quartz cell was cooled to about  $100$  °C, about 2 mL of water was added and heated at about  $100 - 150$  °C for 2 h, and volume was made in 10 mL flask.


## RESULT AND DISCUSSION

The complex of Nickel and Cobalt with dimethylglyoxime ligands (DMG) was formed in an ammonium buffer media.

$\text{Ni}(\text{H}_2\text{O})_6^{2+} + 2 \text{H}_2\text{DMG} \rightarrow \text{Ni}(\text{HDMG})_2 + 2 \text{H}^+ + 6 \text{H}_2\text{O}$ .  
Accumulation take place at control potential, Ni or Co complex was adsorbed at the working electrode surface [4]



In a stripping step, potential was scanned in negative direction (cathodic stripping) so that electrochemical reduction of adsorbed complex of Ni or Co was take place, and formed peak reduction voltammogram.



Current peak (Ip) of cathodic stripping Ni/Co-DNG complex was used as the based of Ni/Co concentration calculation.

Fig 1 shows the dependence of current peak (Ip) of Co/Ni -DMG complex on the ligands concentration. Both current peak (Ip) of Ni and Co complex maximum on the  $5 \cdot 10^{-4}$  M DMG concentration. Therefore this concentration was selected as the optimum concentration of DMG.



Accuracy and reliability of a simultaneously two elements analysis in voltammetry depend on that two peak potential (Ep) separation. According to Polarographic analyzer manual [6], if the peak potential separation more than 100 mV the two elements has no interference each other. Table 1 shows the dependence of Ep on the buffer concentration, from 0.01 M to 1.0 M. The separation between two peaks Ni and Co less than 100 mV only shows at using 0.01 M buffer concentration. While Table 2 shows the effect of pH supporting electrolyte in the range of pH 6 to 11, there are only at pH 9 and 10 give more than 100 mV separation peak. There for using buffer concentration upper than 0.01 M and pH 9 or 10 was safe for simultaneous Ni and Co determination. Fig 2. shows the effect of buffer pH as a supporting electrolyte on current peak (Ip) Ni and Co. The optimum Ip take place on the pH 8-9, and became drop at pH 10.

**Table 1.** Peak potential of Ni and Co in the buffer concentration.

Buffer concentration (M)	Peak potential Ep (V vs Ag/AgCl)		Peak potential separation (mV)
	Ni	Co	
0.01	-0.992	-1.090	98
0.1	-1.022	-1.146	124
0.5	-1.024	-1.166	142
1.0	-1.022	-1.178	156

**Table 2.** Peak potential of Ni and Co in pH buffer variation

pH buffer	Peak potential, Ep (V vs Ag/AgCl)		Peak potential separation (mV)
	Ni	Co	
6	-0.984	-0.994	10
7	-0.984	-0.996	12
8	-0.984	-1.016	32
9	-0.982	-1.092	110
10	-1.032	-1.178	146
11	-1.050	-1.106	56



**Fig 2.** The effect of pH buffer on current peak (Ip) Ni and Co. DMG concentration  $5 \times 10^{-4}$  M, deposition potential  $-0.7$  V, deposition time 120 s, Ni concentration = Co concentration = 20 ng/mL.**Fig 3.** The effect of deposition potential on current peak (Ip) Ni and Co. DMG concentration  $5 \times 10^{-4}$  M, deposition time 120 detik, Ni concentration = Co concentration = 20 ng/mL.

So, the buffer concentration 0.1 – 1 M at pH 9 was selected as the optimum condition of supporting electrolyte for Ni and Co simultaneous determination.

The effect of potential deposition on current peak of Ni and Co was studied from  $-0.5$  V to  $-0.9$  V. In Fig 3 shows that Ip Co is relative constant but for Ip Ni tend to decrease after  $-0.7$  V deposition potential. Therefore,  $-0.7$  V was selected as the optimum deposition potential. Fig 4. shows the effect of deposition times from 60 to 600 s, on the current peak (Ip) both of Ni and Co. the maximum peak current was observed at 180 s for Ni or Co. Therefore 180 s was selected as the optimum deposition time.

By using  $-0.7$  V deposition potential will give opportunity for the other ions like Cu, Cd or Pb were deposited on the surface of working electrode. Because of that Ni and Co determination using the negative scanning (cathodic stripping), Cu, Cd and Pb ion have no interference.

From the result above, we chose the following as the optimum experimental condition: supporting electrolyte is (0.1 – 1) M buffer ammonium tartaric pH  $9.5 \times 10^{-4}$  M DMG,  $-0.7$  V deposition potential and 180 s

**Fig 4.** The effect of deposition time on current peak (Ip) Ni and Co. DMG concentration  $5 \times 10^{-4}$  M, deposition potential  $-0.7$  V, Ni concentration = Co concentration = 20 ng/mL.**Fig 5.** The relationship between current peak (Ip) and Ni and Co concentration in tartaric buffer pH 9, DMG concentration  $5 \times 10^{-4}$  M, deposition potential  $-0.7$ , deposition time 180 s.

**Table 3.** Ni and Co in standard materials (n = 7 – 10)

Material standard	Ni ( $\mu\text{g/g}$ )		Co ( $\mu\text{g/g}$ )	
	Our result	Certificate	Our result	Certificate
Apple leave SRM 1515	0.90 $\pm$ 0.1	0.90 $\pm$ 0.1	nd	nd
Bovine liver SRM 1577B	0.219 $\pm$ 0.008	Nd	0.143 $\pm$ 0.005	nd
Oyster Tissue SRM 1566b	1.047 $\pm$ 0.029	1.04 $\pm$ 0.09	0.361 $\pm$ 0.019	0.371 $\pm$ 0.009
Sargasso CRM no 9	Nd	nd	0.121 $\pm$ 0.005	0.120 $\pm$ 0.01
Copepod RMMAA-1	1.89 $\pm$ 0.17	1.59 $\pm$ 0.2	0.125 $\pm$ 0.015	0.12 $\pm$ 0.01

nd = no data

**Table 4.** Recovery and RSD of Ni and Co in reference materials

Material standard	Recovery (%)		RSD(%)	
	Ni	Co	Ni	Co
Apple leave SRM 1515	100	nd	11	nd
Bovine liver SRM 1577B	nd	nd	3.65	3.50
Oyster Tissue SRM 1566b	100.67	97.30	2.77	5.37
Sargasso CRM no 9	nd	100.8	nd	4.13
CopepodRMMAA-1	103.15	104.17	10.37	6.40

**Table 5.** "t" value calculation and table [8] of the determination result

Material standard	"t" calculation		"t" table
	Ni	Co	
Apple leave SRM 1515	0	nd	2.31
Bovine liver SRM 1577B	nd	nd	-
Oyster Tissue SRM 1566b	0.11	0.182	2.57
Sargasso CRM no 9	nd	0.037	2.57
Copepod RMMAA-1	2.3	0.127	2.31

**Table 6.** Ni and Co contained in the local vegetables

Sample	Ni ( $\mu\text{g/g}$ )	Co ( $\mu\text{g/g}$ )
Cassava leave	0.3 $\pm$ 0.1	nd
Spinach	1.3 $\pm$ 0.2	0.1 $\pm$ 0.0
Mustard greens	nd	nd
Leafy vegetables	nd	nd
carbbage	nd	nd

deposition time. From this optimum experimental condition, a linier relationship between the current peak and Ni or Co concentration was obtained in the range 5 – 30 ng/mL (Fig 5). Linier regression equation for standard curve in this range for Ni is  $Y = 477.9 X + 3548$  with  $r = 0.9998$  and for Co is  $Y = 388.3 X + 228.7$  with  $r = 0.9991$ . Limit detection of Ni and Co is 0.4 ng/ml, obtained by that statistical calculation of standard curve [7]. An application of this method in the determination of Ni and Co in biological material standard showed in Table 3, 4 and 5. Shows that recovery of this method is in the range of 97.30% - 104.17% with relative standard deviation (RSD) is in the range of 2.77% – 10.37%. From the 't' test [7], shows that "t" calculation was smaller than "t" table value.

Therefore, it can be concluded that concentration determination result and the standard concentration value have no significant difference. This result demonstrated that digestion sample have been done perfectly.

## CONCLUSION

The optimum experimental condition: supporting electrolyte was (0.1 – 1 ) M buffer ammonium tartaric pH 9,  $5 \times 10^{-4}$  M DMG, -0.7 V deposition potential and 180 s deposition time. A linier relationship between the current peak and Ni or Co concentration was obtained in the range 5 – 30 ng/mL. Limit of detection of this method 0.4 ng/mL for Ni and Co. The recovery is about 97.30 -104.17% with relative standard deviation (RSD) is about 2.77 -10.37%. Digestion sample using the mixture of  $\text{HNO}_3$ , perchloric acid,  $\text{H}_2\text{O}_2$  were heated to dryness was suitable to the Ni and Co in the biological material determination. There are only cassava leave and spinach that contain Ni and Co in all local vegetables that studied

## REFERENCES

1. Bolt A.W., 1995, *Curr. Separations.*, 14, 1
2. Wang Y, 1985, Stripping analysis. Principles. Instrumentation and application
3. Neiman E.Y and Dracheva L.V, 1990, *Z. Anal. Khim.*, 45, 2
4. Vos L., 1986, *Anal. Chem. Acta.*, 184
5. Thunus L., 1990, *Anal. Lett.* 23, 5
6. Anonim, 1988, *Polarographic Analyzer Instruction Manual, Model 384B*, Princeton Applied Research Corporation
7. Miller J.C and Miller J.N., 1998, *Statistics for Analytical Chemistry*, Ellis Horwood Limited, London